Puławy, 2012.07.04

ZHZ-066/490/12

## OPINIA

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AOAC Research Institute. Certificate of Performance Tested <sup>SM</sup> Status. Certificate No. 011102. Delvotest SP-NT Visual

Bilandžić N., Solomun Kolanović B., Varenina I., Scortichini G., Annunziata L., Brstilo M., Rudan N.: Veterinary drug residues determination in raw milk in Croatia. *Food Control*, 2011, 22, 1941-1948

da S. Costa A., Lobato V.: Avalicao da presenca de residuos de antimicrobianos em leite e bebida lactea UHT por teste de inibicao microbiana comercial. *Rev I Latic „Candida Tostes”*, 2009, 64, 72-76

Le Breton M.H., Savoy-Perroud M.C., Diserens J.M.: Validation and comparison of the Copan Milk Test and Delvotest SP-NT for the detection of antimicrobials in milk. *Anal Chim Acta*, 2007, 586, 280-283

Manafi M., Hesari J., Abbas Rafat S.: Monitoring of antibiotic residue in raw and pasteurised milk in East Azerbaijan of Iran by Delvotest method. *J Food Res*, 2011, 20/3, 125-131

Navratilova P.: Screening methods used for the detection of veterinary residues in raw cow milk – a review. *Czech J Food Sci*, 2008, 26, 393-401

Nikolić N., Mirecki S., Blagojević M.: Presence of inhibitory substances in raw milk in the area of Montenegro. *Mljekarstvo*, 2011, 61, 182-187

NSW Food Authority. General Circular 03/2006. Antibiotic detection method approved for use.

O'Donnell S.: The efficacy of antibiotic residue screening test for the detection of natural antimicrobials in milk. Academic Dissertation. University of Connecticut, 2011

Salomskiene J., Jonkuvienė D., Paserpskiene M.: Comparison of Delvo-x-press and other methods for determining inhibitors in milk. *Vet Med. Zool*, 2009, 45, 73-80

Stead S.L., Ashwin H., Richmond S.F., Sharman M., Langev Stark P.C., Keely B.J.: Evaluation and validation according to international standards the Delvotest SP-NT screening assay for antimicrobial drugs in milk. *Int Dairy J*, 2008, 18, 3-11

Syrovkova Goffova Z., Kožarova I., Mate D., Marcincak S., Gondova Z, Sopkova D.: Comparison of detection sensitivity of five microbial inhibition tests for the screening of aminoglycoside residues in fortified milk. Czech J Food Sci, 2012, 4, 314-320

Wang J.: Analysis of antibiotics in milk and its products. Chapter 35 s. 887-905. CRC Press, 2010. In: Safety analysis of foods of animal origin.

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